Instructions for use of the Bowie & Dick Type Autoclave Lean Pack

Conform EN-ISO 11140-1 (Type 2), EN-ISO 11140-4

UK

This B&D Type Test should be used in the same manner as the standard B&D test pack (EN285). Prior to performing the test, run a 134°C cycle with the chamber empty. Immediately after this cycle place the B&D Type Test in the chamber as closely as possible to the position that during validation is identified to be the cold spot. Undertake a steam penetration test (134°C, max. 3.5 min) in an otherwise empty chamber. Upon completion of the cycle remove the B&D Type Test from the sterilizer, open it and take out the indicator sheet. The sheet must have changed from light blue to pink. If any difference in this colour is noted, repeat the test. If also this second test fails it should be recorded as fail on the testsheet. Record the results on the test sheet with other relevant data. In case of fail the sterilizer is not functioning properly and needs maintenance, adjustment and/or subsequent validation. The indicator sheet should be included in the sterilization file. Always place the Bowie & Dick Lean Test Pack horizontally in the sterilizer and on its back. If you place the pack in a different position in the sterilizer, there is a chance that the test could give a false fail.

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