



# For Radiodermatitis

Soothing and healing with Flamigel<sup>4,5</sup>

# Clinically proven treatment for radiodermatitis<sup>1,2</sup>

Radiodermatitis occurs in 95% of patients after radiotherapy.<sup>1</sup>

## The skin can:

- Become red
- Feel painful
- Itch
- Have a warm or burning sensation

## Patient benefits:

- Flamigel® delays and reduces occurrence of radiotherapy-induced moist desquamation<sup>1,2</sup>
- Cooling effect eases the pain and irritation<sup>2,3</sup>
- Easy to use
- Easy to remove
- No bandage required
- Does not stain clothes
- Aids recovery from the delayed skin reaction to radiotherapy<sup>1,2</sup>

Radiodermatitis can occur 1 to 2 weeks after the radiation treatment.



Radiodermatitis can cause the dose of radiation to be reduced or for the treatment to be stopped completely.

Cools and soothes the pain.<sup>5</sup>

Aids the skin's natural recovery.<sup>2,3</sup>



**LEVEL 1:**  
Light erythema,  
dry desquamation.




**LEVEL 2:**  
Moderate to severe erythema,  
irregular moist desquamation,  
limited to the folds of the skin



**LEVEL 3:**  
Emerging moist desquamation,  
not limited to the folds of the skin.

## Recommended for moist and dry wounds




### IN MOIST WOUNDS

Hydrocolloids absorb exudate

Ideal wound moisture

Accelerates wound healing



### IN DRY WOUNDS

Hydrogel rehydrates

Ideal wound moisture

Accelerates wound healing

## Efficacy of a hydroactive colloid gel versus historical controls for the prevention of radiotherapy-induced moist desquamation in breast cancer patients.

Censabella S. et al. Eur J Oncol Nursing 2017; 29: 1-7

### OBJECTIVE:

To reduce and delay radiation-induced moist desquamation in breast cancer patients.

### METHOD:




All patients were administered equal doses of radiation.

### FINAL STUDY:

Incidence of radiotherapy-induced moist desquamation.

Time until the occurrence of moist desquamation.



### RESULTS:

Flamigel® brings a significant improvement, both preventively and curatively, in the occurrence of radiotherapy-induced moist desquamation compared to dextrapanthenol.

### CONCLUSION :

Using Flamigel® from the start of the radiotherapy treatment compared to dextrapanthenol results in:

- later occurrence of moist desquamation,
- less incidence of moist desquamation,
- Use of 1 product during radiotherapy is much easier for patients than the use of 2 different products.

# Treatment of radiation damage with Flamigel® and Flaminal® Hydro

Henny Hooiman, Nurse, Careyn, Maasland, Netherlands

## BACKGROUND:

**84-year-old woman with a carcinoma in the mouth. Radiation damage to the neck following 33 treatments.**

Patient was treated with Silver sulfadiazine (applied thickly on the neck) and Cutisorb<sup>®</sup> plus a fixing bandage. The wound area worsened and application of the Silver sulfadiazine was very painful. A new wound care approach began on 3 April after discussion with the radiologist at the Erasmus MC.



### DAY 1: INTAKE

Wound was painful. The exudate had no odour. The colour of the exudate was yellow and green.

- T: 100% red
- I: Painful, swollen, pus
- M: Yellow/green and moist
- E: Scabs

**Administration:** 1x per day change of bandage and rinse. Flamigel<sup>®</sup> 0.5 cm cover with Mepitel<sup>®</sup> One and absorbent bandage.



### DAY 24: INSPECTION BY WOUND CONSULTANT

Wound healing successfully. Wound area was much smaller.

- T: 80% red
- I: No and less pain
- M: Very moist
- E: Scab formation, dry

**Administration:** continue current regime.




### DAY 31: INSPECTION BY WOUND CONSULTANT

Moisture in the wound reducing.

- T: 100% red
- I: None, no pain, itchiness
- M: Almost dry
- E: Dry

**Administration:** White islands treated with Flaminal<sup>®</sup> Hydro and remaining skin with Flamigel<sup>®</sup>



### DAY 49: INSPECTION BY WOUND CONSULTANT

Wound closed.

**CONCLUSION:** Recovery of the wound caused by radiation damage was successful with Flamigel<sup>®</sup> and Flaminal<sup>®</sup>. The patient experienced no pain when the dressing was changed.

**INDICATIONS:** Flamigel<sup>®</sup> is indicated as a treatment of superficial burns, including those caused by radiotherapy, as well as for the treatment of minor wounds and superficial open wounds. Flaminal<sup>®</sup> Hydro is indicated to be used on slightly to moderately exuding wounds.

#### References:

- Censabella S. et al. Efficacy of a hydroactive colloid gel versus historical controls for the prevention of radiotherapy-induced moist desquamation in breast cancer patients. Eur J Oncol Nurs. 2017 Aug; 29:1-7.
- Censabella S. et al. Retrospective study of radiotherapy-induced skin reactions in breast cancer patients; reduced incidence of moist desquamation with a hydroactive colloid gel versus dexpanthenol. Eur J Oncol Nurs. 2014 Oct; 18(5):499-504.
- Van der Plas DL. et al. Treatment of recalcitrant wounds of diverse etiology with a new hydroactive gel. Wounds, 2009 Sep; 21(9):243-8.
- Thomas, S. Wound management and dressings. The pharmaceutical press London. 19990 Ch.2: functions of wound dressing, p9-19.
- Data on file. Tests with regards to absorbency and cooling. 2019

© Registered Trademark. Flamigel<sup>®</sup>, Flaminal<sup>®</sup> Hydro are the registered trademarks of Flen Health SA.

REF	Description	UOM
1005	Flamigel 50g	Each