



## OWEAR® SPLASH RESISTANT GOWNS



Provides arm and body protection, ensuring coverage at all times. AAMI Level 1 barrier protection.

### PRODUCT FEATURES

- o One size fits most
- o Waterproof
- o Two back ties to fit comfortably
- o Two neck ties to prevent hole loop slipping

### SUITABLE FOR

- o Hospital and medical
- o Emergency services
- o Mortuary services
- o Aged care



Blue



White



Yellow

REORDER	SIZE	FEATURE	COLOUR	PCS/UNIT	UNITS/CTN
GN001	One Size	AAMI Level 1	Blue	50	1
GN002	One Size	AAMI Level 1	White	50	1
GN003	One Size	AAMI Level 1	Yellow	50	1

All product information is correct at date of publication

