

## SENTURIAN® BASIC DRESSING PACKS



Convenient and easy to use, Sentry Medical dressing packs provide an all in one solution for basic wound care needs. Requires no preparation time and are sterile for cleanliness.

### PRODUCT FEATURES

- o All-in-one pack
- o Easy to use

### SUITABLE FOR

- o Hospital and medical
- o First aid
- o Minor wounds
- o Cuts and grazes

### COMPONENTS

- o PE/Non-woven Sterile Field, Blue (50 x 50cm)
- o White Standard Dressing Tray Type 3 (140mm x 130 mm x 18mm) - compartments: small 30ml, large 90ml
- o White Graded Dressing Tray Type 4 (140 x 130 x 18mm) - compartments: small 20ml, large 80ml
- o Dressing Towel (40cm x 40cm)
- o Dressing Towel Polylined (50cm x 50 cm)
- o Cotton balls
- o Non Woven Balls
- o Gauze Swabs - 8ply (7.5cm x 7.5cm)
- o Non Woven Swabs - 4ply (7.5cm x 7.5cm)
- o Disposable Plastic Forceps: 2 x Blue, 1 x yellow (11cm)

All product information is correct at date of publication



## SENTURIAN® BASIC DRESSING PACKS

REORDER >	WDP004	WDP005	WDP006	WDP007	WDP009	WDP010	WDP011	WDP013	WDP014	WDP021
Opening Type	Tear	Tear	Tear	Peel	Peel w/ Tuck-in Apex	Tear	Tear	Tear	Peel	Peel
PE/Non-woven Sterile Field	•	•	•	•	•	•	•	•	•	•
White Standard Dressing Tray	•	•					•			
White Graded Dressing Tray										
Dressing Towel	•	•	•	•	•	•	•	Polylined	Polylined	Polylined
Cotton balls	5pcs									
Non Woven Balls		5pcs								
Gauze Swabs	2pcs			6pcs			5pcs		6pcs	
Non Woven Swabs		2pcs	5pcs		6pcs	6pcs		6pcs		6pcs
Disposable Plastic Forceps	•	•	•	•	•	•	•	•	•	•



M E D I C A L

1300 995 999 | sales@sentrymedical.com.au  
www.sentrymedical.com.au