

ThreadX[®]

THREADX[®] X-RAY DETECTABLE NON-WOVEN BALLS



ThreadX[®] Non-Woven Balls are constructed from a polyester and rayon blend, with a black x-ray detectable thread heat bonded to the material.

PRODUCT FEATURES

- o X-ray detectable
- o Black thread for ease of visibility
- o Highly absorbent
- o Low linting

SUITABLE FOR

- o Open surgical procedures
- o Absorption of exudate



REORDER	PCS	PCS/UNIT	UNITS/CTN
XNB020	5pcs	600	1

All product information is correct at date of publication



1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

BR10081_V1_11_18