

# ThreadX<sup>®</sup>

## THREADX<sup>®</sup> X-RAY DETECTABLE GAUZE




ThreadX<sup>®</sup> Optimum Gauze Swabs are 100% cotton, with various weave count options for preferred absorbency. A black x-ray detectable thread is heat bonded to the gauze.

### PRODUCT FEATURES

- o French folded to reduce fraying
- o Low linting
- o X-ray detectable
- o Highly absorbent
- o Optimum and Standard weave count variants

### SUITABLE FOR

- o All open surgical procedures
- o Operating room



Sterile Packaging

REORDER	PLY	SIZE	WEAVE COUNT	PCS	STERILITY	PCS/UNIT	UNITS/CTN
XRD020	8ply	10cm x 10cm	Optimum	5pcs	Non Sterile	600	1
XRD3252	8ply	10cm x 10cm	Standard	5pcs	Sterile	250	1
XRD3253	8ply	10cm x 10cm	Optimum	5pcs	Sterile	250	1

All product information is correct at date of publication

