

# ThreadX<sup>®</sup>

## THREADX<sup>®</sup> X-RAY PACKING GAUZE



ThreadX<sup>®</sup> Packing Gauze is constructed from 100% cotton. Designed for open surgical procedures, an x-ray detectable thread is heat bonded to the gauze.

### PRODUCT FEATURES

- o X-ray detectable thread
- o Absorbs fluid and protects organ and tissue injury
- o Helps control bleeding
- o Low linting

### SUITABLE FOR

- o Open surgical procedures
- o Deep wounds
- o Undermining and tunnels
- o Heavy exudate procedures (XRV002, XRV003)



REORDER	PLY	SIZE	THREAD COLOUR	PCS/UNIT	UNITS/CTN
XRP001	4ply	1.25cm x 1m	Blue	100	1
XRP006	4ply	2.5cm x 1m	Blue	100	1
XRP007	4ply	2.5cm x 2m	Blue	100	1
XRP011	4ply	5cm x 1m	Blue	100	1
XRP012	4ply	5cm x 2m	Blue	100	1
XRP016	4ply	10cm x 1m	Blue	100	1
XRP017	4ply	10cm x 2m	Blue	100	1
XRV002 <sup>&lt;&lt;</sup>	8ply	10cm x 2.5m	Blue	100	1
XRV003	8ply	10cm x 3m	Black	100	1

<< Indent Item

All product information is correct at date of publication

