

Eye Pad

Asguard Eye Pads are constructed with an absorbent inner and covered with a non-woven outer layer with sealed edges designed to protect the eye or absorb secretions.



Shaped to Fit

Oval shaped to specifically fit over the eye socket



Sealed Edges

Non-woven outer layer with sealed edges prevents linting



Soft and Comfortable

Non-abrasive outer layer soft on delicate skin



Suitable for:

- Minor eye injuries
- Post-surgical dressing
- First aid kits

REF	Size	Sterility	Pcs/Unit	Units/Ctn
EP001	6cm x 8cm	Sterile	50	24
EP002	6cm x 8cm	Non-sterile	2000	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

