

Transparent Film Dressing Roll

The Asguard Film Dressing Roll is a transparent adhesive film with moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. Asguard Film Dressings conform to the skin and allow for easy visibility to wound healing.



Water Resistant

Impermeable to help protect site against infection and allows showering with the dressing in-situ



Cut to Length

Measured grid backing allows cut to correct size as required



Contours to Body

Flexible and comfortable layer adheres to the shape of application site



Suitable for:

- Wounds with little or no exudate
- Protection and securement of IV lines
- Dressing reinforcement

REF	Size	Pcs/Unit	Units/Ctn
ISB010	5cm x 10m	1	24
ISB011	10cm x 10m	1	12

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

