

Conforming Bandage

Bodichek Conforming Bandages are specifically designed to adapt to the shape of the body area applied to, ensuring a secure fit without restricting movement or circulation.

Lightweight and breathable, conforming bandages are comfortable to wear for extended periods and are suitable for a wide range of medical applications, including post-surgical care, sprain management and protecting wound dressings.



Breathable

Porous non-woven material provides high comfort and air circulation



Cut to Length

Non-fraying edges when cut to length. Can be affixed with tape or clips.



Conforms to Body Contours

Allows flexibility without restriction of movement



Suitable for:

- Dressing retention
- Light compression and support

REF	Size	Sterility	Pcs/Unit	Units/Ctn
CFB001	2.5cm x 1.5m (Unstretched)	Non-Sterile	12	120
CFB002	5cm x 1.5m (Unstretched)	Non-Sterile	12	60
CFB003	7.5cm x 1.5m (Unstretched)	Non-Sterile	12	40
CFB004	10cm x 1.5m (Unstretched)	Non-Sterile	12	30
CFB005	15cm x 1.5m (Unstretched)	Non-Sterile	12	20

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

