

Tourniquet

The Bodichek Tourniquet consists of a strap and clip buckle mechanism to prevent pinching, catching or jamming. These tourniquets feature a button that gradually releases pressure without opening, allowing the band to be loosened and re-tightened with one-handed operation. Designed to be used multiple times on a single patient, to help prevent cross infections.



Easy to Use

Pressure release button allows for one-handed loosening and tightening as well as easy application and removal



Comfortable Band

Soft and adjustable band to suit every patient



Single Patient Use

Helps prevent risk of patient infection



10180001 Dispenser



10180001 Bag



TNQ001



10180001

- Suitable for:**
- Blood sample collection
 - Cannulation procedure
 - First aid use to stop bleeding

REF	Size	Band Thickness	Colour	Feature	Dispenser Box	Sterility	Pcs/Unit	Units/Ctn
TNQ001	2.5cm x 40cm	1.95mm	Yellow band	Clip buckle mechanism	No	Non-Sterile	100	5
10180001	2.5cm x 40cm	2.05mm	Blue band	Clip buckle mechanism	Yes	Non-Sterile	50	6

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
 1300 995 999 | sales@sentrymedical.com.au
 www.sentrymedical.com.au

