

Tubular Support Bandage 1 Metre

Bodigrip Tubular Support Bandage is an elasticated bandage designed to provide support for soft tissue injuries such as swelling, bruising and weak joints. The seamless knitted design provides even support without restricting movement and eliminates the need for tapes or pins.



Conforms to Body Contours

Designed with light elasticity for a comfortable and conformable fit on various body parts



Reusable

Can be washed and reused, while maintaining integrity and support during extended use



Extensive Range

Available in a range of sizes to suit all limbs and joints



- Suitable for:**
- Support for soft tissue injuries such as strains and sprains
 - Light to moderate compression to alleviate swelling
 - Securing dressings

REF	Size	Bandage Size	Colour	Application	Pcs/Unit	Units/Ctn
13023161	B	6.5cm x 1m	Flesh	Small wrist/elbow	12	6
13023162	C	6.75cm x 1m	Flesh	Medium wrist/elbow, small ankle	12	6
13023163	D	7.5cm x 1m	Flesh	Large wrist/elbow, medium ankle, small knee	12	6
13023164	E	8.75cm x 1m	Flesh	Large ankle, medium knee, small thigh	12	6
13023165	F	10cm x 1m	Flesh	Large knee, medium thigh	12	6
13023166	G	12cm x 1m	Flesh	Large thigh	12	6

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

