

IV Starter Packs Overview

Convenient and easy to use, the Sage IV Starter Pack provides an all-in-one solution with essential items to initiate intravenous therapy swiftly and effectively.

The IV Starter Pack helps maintain a sterile environment, ensuring that healthcare professionals can perform the procedure safely and hygienically, reducing the risk of infections.



All-in-one Pack

Comprehensive components help streamline medical procedures



Aids in Infection Control

Provides an on-demand sterile set-up to help minimise risk of infection



Clinical and Time Efficient


Access to vital components for IV therapy and easy setup



IV Starter Pack #9

Suitable for:

- Outpatient clinics
- Surgical procedures
- Emergency rooms and trauma care




Comes in a convenient desktop dispenser of 45 pieces.

REF	Description	Sterility	Pcs/Unit	Units/Ctn
11005042	IV Starter Pack #9	Sterile	180	1
11005044	IV Starter Pack #11	Sterile	180	1
11005045	IV Starter Pack #12*	Sterile	180	1
11005046	IV Starter Pack #13*	Sterile	180	1

*Minimum order quantities apply.

Component details

IV Starter Packs Overview



ISB301

(IV Starter Pack #12)




ISB302

(IV Starter Pack #11)



ISB303

(IV Starter Pack #9)



ISB304

(IV Starter Pack #13)

REF	Description	Components	Sterility	Pcs/Unit	Units/Ctn
11005042	IV Starter Pack #9	<ul style="list-style-type: none"> 1 x Blue Polylined Towel - 42cm x 40cm 1 x IV Starter Pack Tray 1 x Swab Stick, 2% CHG, 70% Isopropyl Alcohol 3 x Non-Woven Swabs - 7.5cm x 7.5cm, 4ply 1 x Asguard IV Dressing Reinforced 7cm x 8.5cm (ISB303) 2 x Patient Label - Green Relocated IV Site 1 x White Plastic Field 47.5cm x 47.5cm 	Sterile	180	1
11005044	IV Starter Pack #11	<ul style="list-style-type: none"> 1 x Blue Polylined Towel - 42cm x 40cm 1 x IV Starter Pack Tray 1 x Swab Stick, 2% CHG, 70% Isopropyl Alcohol 3 x Non-Woven Swabs - 7.5cm x 7.5cm, 4ply 1 x Asguard IV Dressing Reinforced 6.5cm x 7cm (ISB302) 2 x Patient Label - Green Relocated IV Site 1 x White Plastic Field 47.5cm x 47.5cm 	Sterile	180	1
11005045	IV Starter Pack #12	<ul style="list-style-type: none"> 1 x Blue Polylined Towel - 42cm x 40cm 1 x IV Starter Pack Tray 1 x Swab Stick, 2% CHG, 70% Isopropyl Alcohol 3 x Non-Woven Swabs - 7.5cm x 7.5cm, 4ply 1 x Asguard IV Dressing Reinforced 5cm x 5.7cm (ISB301) 2 x Patient Label - Green Relocated IV Site 1 x White Plastic Field 47.5cm x 47.5cm 	Sterile	180	1
11005046	IV Starter Pack #13	<ul style="list-style-type: none"> 1 x Blue Polylined Towel - 42cm x 40cm 1 x IV Starter Pack Tray 1 x Swab Stick, 2% CHG, 70% Isopropyl Alcohol 3 x Non-Woven Swabs - 7.5cm x 7.5cm, 4ply 1 x Asguard IV Dressing Reinforced 7cm x 9cm (ISB304) 2 x Patient Label - Green Relocated IV Site 1 x White Plastic Field 47.5cm x 47.5cm 	Sterile	180	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

