

Bed Draw Sheet

Owear Bed Draw Sheets are impervious, and treated with an anti-static solution for ease of use and maximum patient comfort.



Infection Control

Provides a protected surface in clinical settings



Convenient

Easy to use with simple application



Anti-Static

Enhancing patient comfort



Suitable for: - Examination beds

REF	Size	Colour	Sterility	Pcs/Unit	Units/Ctn
DS001	91.5cm x 152cm	Blue	Non-sterile	300	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

