

Face Shield

A protective face shield to help protect the eyes, nose, and mouth from sprays and splatters.



Anti Fog Coating

Clear vision with dual-sided anti-fog treatment and protective film



Universal Fit

Elastic band accommodates all head sizes



Comfortable

Soft forehead cushion ensures comfort during prolonged wear



- Suitable for:**
- Sterilisation department
 - Emergency rooms and isolation wards
 - Splash and dust protection

REF	Size	Sterility	Pcs/Unit	Units/Ctn
10025515	25cm x 29cm	Non-Sterile	200	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

