

Impervious Gowns

Owear Impervious Gowns provide effective protection against fluids with both the front and back of the gown being impervious. Gowns meet AAMI Level 2 liquid barrier performance levels, ensuring a reliable level of fluid resistance suitable for low to moderate risk scenarios.



AAMI Level 2

Gowns provide low barrier against fluid exposure and contaminants



Chemotherapy Tested *

Provides protection when handling and administering chemotherapy **




Neck and Waist Ties


Elastic cuff with neck and back ties to easily don and doff ***



Blue



Yellow



Purple
(Oncology)

Suitable for:

- Hospital and medical
- Laboratory
- Low risk patient isolation scenarios
- Light exposure to blood or body fluids
- Pathology and ICU

REF	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN004P	Regular	AAMI Level 2	Purple (Oncology)	100	1
GN005P	Regular	AAMI Level 2	Blue	50	1
GN006P	Regular	AAMI Level 2	Yellow	50	1

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

*** GN004P - velcro neck tie.

Component details

REORDER	AAMI LEVEL 1			AAMI LEVEL 2							AAMI LEVEL 3		Thumb Loop		
	GN001	GN002	GN003	GN001P	GN036P	GN003P	GN056P	GN0041P	GN004P	GN005P	GN006P	GN009P	GN334	GN336	GN017P
Size	Regular	Regular	Regular	Regular	X-Large	Regular	X-Large	Regular	Regular	Regular	Regular	Regular	Universal	X-Large	Large
Colour	Blue	White	Yellow	Blue	Blue	Yellow	Yellow	Purple	Purple	Blue	Yellow	White	Blue	Blue	Blue
Sterility	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile
Pcs/Unit	50	50	50	50	50	50	50	50	100	50	50	50	50	50	120
Unit/Ctn	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Back Ties	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Neck Ties	•	•	•							•	•	•			
Velcro Neck				•	•	•	•	•	•				•	•	
Halter Neck															
Polybreath												•			•
Sof-Cuff				•	•	•	•	•	•				•	•	
Elastic Cuff	•	•	•							•	•	•			
Thumb Loop															•
Short Sleeve															
Impervious				with gloves	with gloves	with gloves	with gloves	with gloves	•	•	•	•	with gloves	with gloves	
Splash Resistant	•	•	•												
Chemically Tested**				•	•	•	•	•	•	•	•	•			

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.