

Impervious Gowns - Sof-Cuff

Owear Impervious Gowns - Sof-Cuff provide effective protection against fluids with the front and back of the gown being impervious. Owear Impervious Gowns meet the AAMI Level 2 liquid barrier performance levels, ensuring a reliable level of fluid resistance suitable for low to moderate risk scenarios.



AAMI Level 2

Gowns provide low barrier against fluid exposure and contaminants






Chemotherapy Tested *

Provides protection when handling and administering chemotherapy **



Velcro Neck and Waist Ties

Velcro neckline and back ties to easily don and doff



Suitable for:

- Hospital and medical
- Laboratory
- Light exposure to blood or body fluids
- Low risk patient isolation scenarios
- Pathology and ICU

REF	Size	Feature	Colour	Sterility	Pcs/Unit	Units/Ctn
GN001P	Regular	AAMI Level-2	Blue	Non-Sterile	50	1
GN003P	Regular	AAMI Level-2	Yellow	Non-Sterile	50	1
GN0041P	Regular	AAMI Level-2	Purple	Non-Sterile	50	1
GN036P	X-Large	AAMI Level-2	Blue	Non-Sterile	50	1
GN056P	X-Large	AAMI Level-2	Yellow	Non-Sterile	50	1

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

Component details

REORDER	AAMI LEVEL 1			AAMI LEVEL 2							AAMI LEVEL 3		Thumb Loop		
	GN001	GN002	GN003	GN001P	GN036P	GN003P	GN056P	GN0041P	GN004P	GN005P	GN006P	GN009P	GN334	GN336	GN017P
Size	Regular	Regular	Regular	Regular	X-Large	Regular	X-Large	Regular	Regular	Regular	Regular	Regular	Universal	X-Large	Large
Colour	Blue	White	Yellow	Blue	Blue	Yellow	Yellow	Purple	Purple	Blue	Yellow	White	Blue	Blue	Blue
Sterility	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile
Pcs/Unit	50	50	50	50	50	50	50	50	100	50	50	50	50	50	120
Unit/Ctn	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Back Ties	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Neck Ties	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Velcro Neck				•	•	•	•	•	•				•	•	
Halter Neck															
Polybreath											•				•
Sof-Cuff				•	•	•	•	•	•	•	•	•	•	•	•
Elastic Cuff	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Thumb Loop															•
Short Sleeve															
Impervious				with gloves	with gloves	with gloves	with gloves	with gloves	•	•	•	•	with gloves	with gloves	
Splash Resistant	•	•	•												
Chemically Tested**				•	•	•	•	•	•	•	•	•			

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

