


L3 Isolation Gowns Soft-Cuff

Owear Isolation Gowns - Sof-Cuff offer protection in environments with moderate levels of fluid exposure. These non-sterile gowns provide a protective barrier with maximum comfort and breathability for the wearer.

3 **AAMI Level 3**
Gowns provide a moderate barrier against fluid exposure and contaminants

 **Velcro Neck and Waist Ties**
Velcro neckline and back ties to easily don and doff

 **Comfortable**
Soft cuffs allow for added comfort and ease when donning gloves



- Suitable for:**
- Hospital and medical
 - Laboratory
 - Light exposure to blood or body fluids
 - Low risk patient isolation scenarios
 - Pathology and ICU

REF	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN334	Universal	AAMI Level-3	Blue	50	1
GN336	X-Large	AAMI Level-3	Blue	50	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

