

Thumb Loop Impervious Gowns

The Owear Overhead Thumb Loop Impervious Gowns ensure arm and body protection from fluids as well as provide an effective barrier against cross infection. Thumb loop helps keep the gown from riding up the arm.



Impervious

An impervious layer helps provide protection against fluids



Easy to Use

Overhead neck closure allows for quick application and easy removal



Everyday Use

Various suitabilities allow for everyday use



Suitable for:

- Hospital and medical
- Laboratory

| REF | Size | Feature | Colour | Pcs/Unit | Units/Ctn |
|--------|-------|----------------------------|--------|----------|-----------|
| GN017P | Large | Overhead Gown - Thumb Loop | Blue | 120 | 1 |

Component details

| REORDER | AAMI LEVEL 1 | | | AAMI LEVEL 2 | | | | | | | AAMI LEVEL 3 | | Thumb Loop | | |
|----------------------------|--------------|-------------|-------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|--------------|-------------|-------------|-------------|-------------|
| | GN001 | GN002 | GN003 | GN001P | GN036P | GN003P | GN056P | GN0041P | GN004P | GN005P | GN006P | GN009P | GN334 | GN336 | GN017P |
| Size | Regular | Regular | Regular | Regular | X-Large | Regular | X-Large | Regular | Regular | Regular | Regular | Regular | Universal | X-Large | Large |
| Colour | Blue | White | Yellow | Blue | Blue | Yellow | Yellow | Purple | Purple | Blue | Yellow | White | Blue | Blue | Blue |
| Sterility | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile | Non-Sterile |
| Pcs/Unit | 50 | 50 | 50 | 50 | 50 | 50 | 50 | 50 | 100 | 50 | 50 | 50 | 50 | 50 | 120 |
| Unit/Ctn | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Back Ties | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Neck Ties | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Velcro Neck | | | | • | • | • | • | • | • | | | | • | • | |
| Halter Neck | | | | | | | | | | | | | | | |
| Polybreath | | | | | | | | | | | • | | | | • |
| Sof-Cuff | | | | • | • | • | • | • | • | • | • | • | • | • | • |
| Elastic Cuff | • | • | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Thumb Loop | | | | | | | | | | | | | | | • |
| Short Sleeve | | | | | | | | | | | | | | | |
| Impervious | | | | with gloves | with gloves | with gloves | with gloves | with gloves | • | • | • | • | with gloves | with gloves | |
| Splash Resistant | • | • | • | | | | | | | | | | | | |
| Chemically Tested** | | | | • | • | • | • | • | • | • | • | • | | | |

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

