

Polybreath Impervious Gowns

Owear Polybreath Impervious Gowns are constructed from an impervious yet breathable material enhancing the wearer's comfort. These Polybreath Gowns meet the AAMI Level 2 liquid barrier performance levels, ensuring a reliable level of fluid resistance suitable for low to moderate risk scenarios.



AAMI Level 2

Gowns provide low barrier against fluid exposure and contaminants



Chemotherapy Tested *

Provides protection when handling and administering chemotherapy **



Sonically Bonded

Strong sealed seams



Suitable for:

- Hospital and medical
- Light exposure to blood or body fluids
- Low risk patient isolation scenarios
- Pathology and ICU
- Laboratory

REF	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN009P	Regular	AAMI Level 2	White	50	1

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

Component details

REORDER	AAMI LEVEL 1			AAMI LEVEL 2							AAMI LEVEL 3		Thumb Loop		
	GN001	GN002	GN003	GN001P	GN036P	GN003P	GN056P	GN0041P	GN004P	GN005P	GN006P	GN009P	GN334	GN336	GN017P
Size	Regular	Regular	Regular	Regular	X-Large	Regular	X-Large	Regular	Regular	Regular	Regular	Regular	Universal	X-Large	Large
Colour	Blue	White	Yellow	Blue	Blue	Yellow	Yellow	Purple	Purple	Blue	Yellow	White	Blue	Blue	Blue
Sterility	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile
Pcs/Unit	50	50	50	50	50	50	50	50	100	50	50	50	50	50	120
Unit/Ctn	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Back Ties	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Neck Ties	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Velcro Neck				•	•	•	•	•	•				•	•	
Halter Neck															
Polybreath											•				•
Sof-Cuff				•	•	•	•	•	•	•	•	•	•	•	•
Elastic Cuff	•	•	•								•	•			
Thumb Loop															•
Short Sleeve															
Impervious				with gloves	with gloves	with gloves	with gloves	with gloves	•	•	•	•	with gloves	with gloves	
Splash Resistant	•	•	•												
Chemically Tested**				•	•	•	•	•	•	•	•	•			

* Chemotherapy drugs tested against is available on request.

** Suitable for use in preparing or handling chemotherapy drugs in oncology departments.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au

www.sentrymedical.com.au

POLYBAG



CARTON

