

Splash Resistant Gowns

Owear Splash Resistant gowns provide arm and body protection, ensuring coverage at all times.



AAMI Level 1

Gowns provide a minimal barrier against fluid exposure and contaminants



Sonically Bonded

Strong sealed seams




Neck and Waist Ties


Adjustable neckline and back ties to secure gown closed



Blue



White



Yellow

Suitable for:

- Daily patient care
- Little to no risk of fluid exposure
- Basic isolation scenarios

REF	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN001	Regular	AAMI Level 1	Blue	50	1
GN002	Regular	AAMI Level 1	White	50	1
GN003	Regular	AAMI Level 1	Yellow	50	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

