

Owear[®] Headwear Overview

Owear[®] headwear range includes a vast assortment of protective head coverings to satisfy the varying demands of healthcare workers. Made from fluid repellent spun-bound polypropylene, the headwear range maintains airflow to ensure comfort for prolonged periods.



Comfortable

Anti-static and low linting for extended wear



Various Colours

Colour coded headwear help to quickly identify health professional's roles in busy settings



Disposable

Assists in infection control maintenance



Round Caps

Suitable for:

- Hospital and medical
- Emergency rooms
- Operating theatres
- Delivery suites
- ICU



Scarves



Theatre Caps



Surgical Hoods

Component details

| REORDER | FEATURE | COLOUR | PCS/ UNIT | UNITS/ CTN | SWATCHES |
|-----------------------|-----------|-----------|--------------|---------------|--|
| ROUND CAPS | | | | | |
| RC001 | Regular | Blue | 250 | 4 |    Blue Green Pink |
| RC001L | Large | Blue | 250 | 4 | |
| RC002 | Regular | Green | 250 | 4 | |
| RC003 | Regular | Pink | 250 | 4 |    White Yellow Christmas |
| RC004 | Regular | White | 250 | 4 | |
| RC005 | Regular | Yellow | 250 | 4 | |
| RC006X | Regular | Christmas | 250 | 1 |   Red Dark Blue |
| RC007 | Regular | Red | 250 | 4 | |
| RC008 | Regular | Dark Blue | 250 | 4 | |
| RCD001 | Dispenser | | 1 | 5 | |
| SCARVES | | | | | |
| SF001 | Regular | Blue | 250 | 2 |    Blue Green Pink |
| SF002 | Regular | Green | 250 | 2 | |
| SF003 | Regular | Pink | 250 | 2 | |
| THEATRE CAPS | | | | | |
| TC005 | Regular | Blue | 250 | 8 |    Blue Green Pink |
| TC006 | Regular | Green | 250 | 8 | |
| TC007 | Regular | Pink | 250 | 8 | |
| TC008 | Regular | White | 250 | 8 | |
| TC009 | Regular | Dark Blue | 250 | 8 | |
| SURGICAL HOODS | | | | | |
| HD001 | Regular | Blue | 250 | 4 |  Blue |
| HD002 | Small | Blue | 250 | 4 | |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.