

Impervious Gowns AAMI Level 2

Owear Impervious Gowns provide effective protection against fluids with the front and back of the gown being impervious. Owear Impervious Gowns meet the AAMI Level 2 liquid barrier performance levels, ensuring a reliable level of fluid resistance suitable for low to moderate risk scenarios.

2 AAMI Level 2
Gowns provide a moderate barrier against fluid exposure and contaminants

Chemotherapy Tested
Provides protection when handling and administering chemotherapy*

Velcro Neck and Waist Ties
Velcro neckline and back ties to easily don and doff**



- Suitable for:**
- Hospital and medical
 - Laboratory
 - Light exposure to blood or body fluids
 - Low risk patient isolation scenarios
 - Pathology and ICU

REF	Size	Feature	Colour	Sterility	Pcs/Unit	Units/Ctn
GN001P	Regular	Soft-Cuff	Blue	Non-Sterile	50	1
GN003P	Regular	Soft-Cuff	Yellow	Non-Sterile	50	1
GN0041P	Regular	Soft-Cuff	Purple (Oncology)	Non-Sterile	50	1
GN00416P	X-Large	Soft-Cuff	Purple (Oncology)	Non-Sterile	50	1
GN036P	X-Large	Soft-Cuff	Blue	Non-Sterile	50	1
GN056P	X-Large	Soft-Cuff	Yellow	Non-Sterile	50	1
GN004P	Regular	Elastic-Cuff	Purple (Oncology)	Non-Sterile	100	1
GN0046P	X-Large	Elastic-Cuff	Purple (Oncology)	Non-Sterile	50	1
GN005P	Regular	Elastic-Cuff	Blue	Non-Sterile	50	1
GN006P	Regular	Elastic-Cuff	Yellow	Non-Sterile	50	1

* Suitable for use in preparing or handling chemotherapy drugs in oncology departments. Chemotherapy drugs tested against is available on request.

** GN005P, GN006P, GN009P - regular neck ties.

Please see Owear Gowns Overview Flyer for full range details.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

