

## Surgical Mask Earloop IIR

Owear Surgical Masks provide a moderate level of protection during surgical procedures. They are 3ply, single use, and create an important barrier for patients against airborne contaminants generated by clinicians and protects clinicians from fluids or particles from the patient.



### Soft and Comfortable

Carefully engineered and clinically preferred for extended periods of use without discomfort



### Disposable

Aids in infection control and minimises the risk on contamination



### Compliance Assurance

Meets the standards of AS 4381:2015 and ISO 14683:2019 Type IIR (120mmHg)



### Suitable for:

- Hospital and medical
- Emergency rooms and operating theatres
- Areas with moderate risk of fluid splash

REF	Quantity	Sterility	Pcs/Unit	Units/Ctn
10025519	50	Non-Sterile	30	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

