



# Sage<sup>®</sup> Basic Dressing Packs

## Classic & Premium Packs


### Aaxis Pacific Aseptic Opening Technique

Non touch opening method ensures pack can be opened without compromising the sterile field

Aseptic transfer point permits tray to be moved without compromising sterility

Rigid squish pad area allows excess fluid to be removed from swabs and bobs and be retained in tray corners

Environmentally safe – incineration process converts pack to harmless gasses of oxygen and carbon dioxide



Non Touch Opening Method

PRODUCT CODE	PRODUCT DESCRIPTION	SKU (MOQ)	SKU/ Inner	SKU/ CTN
<b>CLASSIC DRESSING PACKS</b>				
11001071	Sage <sup>®</sup> Basic Dressing Pack #71 PP MT Sterile , Mini	140	N/A	140
11001074	Sage <sup>®</sup> Basic Dressing Pack #74 PP MT Sterile , Mini	140	N/A	140
11001076	Sage <sup>®</sup> Basic Dressing Pack #76 PP MT Sterile	140	N/A	140
<b>AAT – PREMIUM DRESSING PACKS</b>				
11001001	Sage <sup>®</sup> OPS <sup>™</sup> Basic Dressing Pack #1 PP, Sterile	220	20	220
11001003	Sage <sup>®</sup> OPS <sup>™</sup> Basic Dressing Pack #3 PP, Sterile	220	20	220
11001008	Sage <sup>®</sup> OPS <sup>™</sup> Basic Dressing Pack #8 PP, Sterile	220	20	220
11001010	Sage <sup>®</sup> Basic Dressing Pack #10 PP, Sterile	140	N/A	140
11001114	Sage <sup>®</sup> OPS <sup>™</sup> Basic Dressing Pack #14 PP, Sterile	220	20	220
11001116	Sage <sup>®</sup> OPS <sup>™</sup> Basic Dressing Pack #16 Inverted PP, Sterile	220	20	220
11001068	Sage <sup>®</sup> Basic Dressing Pack #68 PP, Sterile	160	20	160



# Sage<sup>®</sup> Basic Dressing Packs

## Classic - Mini

### Basic Dressing Pack #71

PRODUCT CODE: 11001071

1 x Small Tray 125 x 105 x 25mm

8 x Cotton Ball 0.3g

5 x Non Woven Swabs 7.5 x 7.5cm 4 Ply

2 x Forceps Green

1 x Forceps Yellow

1 x Towel Blue Guard 40 x 40cm

1 x White Plastic Field 50 x 50cm

### Basic Dressing Pack #76

PRODUCT CODE: 11001076

1 x Aseptic Tray

5 x Non-Woven Bobs .3gm

2 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply

2 x OPS Forceps Green

1 x OPS Forceps Yellow

1 x Towel Blue Guard 40 x 40cm

1 x White Plastic Field 50 x 50cm

### Basic Dressing Pack #74

PRODUCT CODE: 11001074

1 x Small Tray

6 x Gauze Swabs 7.5 x 7.5cm 8-Ply

2 x Forceps Green

1 x OPS Forceps Yellow

1 x Towel Blue Polylined 42 x 42cm

1 x Field White Plastic 47.5 x 47.5cm



# Sage<sup>®</sup> Basic Dressing Packs

## AAT – Premium

### Basic Dressing Pack #1

PRODUCT CODE: 11001001

1 x OPS<sup>™</sup> Premium Tray  
5 x Cotton Balls  
2 x Gauze Swabs 7.5 x 7.5cm 8-Ply  
2 x Forceps Green  
1 x Forceps Yellow  
1 x Towel 2-Ply Paper 40 x 40cm  
1 x Field White Plastic 47.5 x 47.5cm

### Basic Dressing Pack #3

PRODUCT CODE: 11001003

1 x OPS<sup>™</sup> Premium Tray  
5 x Gauze Swabs 7.5 x 7.5cm 8-Ply  
2 x Forceps Green  
1 x OPS Forceps Yellow  
1 x Towel 2-Ply Paper 40 x 40cm  
1 x Field White Plastic 47.5 x 47.5cm

### Basic Dressing Pack #8

PRODUCT CODE: 11001008

1 x OPS<sup>™</sup> Premium Tray  
5 x Non-Woven Bobs .3gm  
2 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply  
2 x OPS OPS Forceps Green  
1 x OPS Forceps Yellow  
1 x Towel 2-Ply Paper 38 x 40cm  
1 x Field White Plastic 47.5 x 47.5cm

### Basic Dressing Pack #10

PRODUCT CODE: 11001010

1 x Small Tray 125 x 105 x 25mm  
5 x Cotton Balls  
2 x Gauze Swabs 7.5 x 7.5cm 8-Ply  
3 x OPS Forceps Green  
1 x Towel 2-Ply Paper 42 x 40cm  
1 x Field White Plastic 47.5 x 47.5cm

### Basic Dressing Pack #14

PRODUCT CODE: 11001014

1 x OPS<sup>™</sup> Premium Tray  
6 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply  
2 x OPS Forceps Green  
1 x OPS Forceps Yellow  
1 x Towel 2-Ply Paper 42 x 40cm  
1 x Field White Plastic 50 x 50cm

### Basic Dressing Pack #16

PRODUCT CODE: 11001116

1 x OPS<sup>™</sup> Premium Tray  
6 x Non-Woven Swabs 7.5cm x 7.5 cm 4-Ply  
2 x Forceps Green  
1 x Forceps Yellow  
1 x Towel Blue Polylined 42 x 40cm  
1 x Field White Plastic 50 cm x 50 cm



# Sage<sup>®</sup> Basic Dressing Packs

## AAT – Premium

### Basic Dressing Pack #44

PRODUCT CODE: 11001144

1 x OPS<sup>™</sup> Premium Tray

6 x Non-Woven Swabs 7.5 x 7.5cm 4 Ply

2 x OPS Forceps Green

1 x OPS Forceps Yellow

1 x Towel 2-Ply Paper 40 x 40cm

1 x Field White Plastic 50 x 50cm

### Basic Dressing Pack #68

PRODUCT CODE: 11001068

1 x OPS<sup>™</sup> Premium Tray

5 x Non Woven Bobs

2 x Non Woven Swabs 7.5 x 7.5cm 4 Ply

2 x OPS Forceps Green

1 x OPS Forceps Yellow

1 x Towel Paper 42 x 40cm

1 x White Plastic Field 50 x 50cm