

## Catheter Packs

Sage Catheter Pack is carefully assembled to ensure convenience, hygiene, and ease of use during catheterisation procedures. They are designed to minimise the risk of infection and provide healthcare professionals with the essential tools needed to perform catheterisation with meticulous care.



### Easy Open Pack

Quick and convenient one-step process



### Puncture Resistant Tray

Firm structure with reinforced edges



### Sterility Guaranteed

Minimising the risk of contamination



11003001



11003009



11003012

### Suitable for:

- Hospitals, clinics, and other healthcare settings where catheterisation is a routine procedure

REF	Description	Sterility	Pcs/Unit	Units/Ctn
<b>Standard Catheter Packs</b>				
11003001	Catheter Pack #1	Sterile	30	1
11003002	Catheter Pack #2	Sterile	60	1
11003009	Catheter Pack #9	Sterile	55	1
<b>OPS Catheter Packs with Purifil Water Syringe</b>				
11003011	OPS™ Catheter Pack #11	Sterile	48	1
11003012	OPS™ Catheter Pack #12	Sterile	48	1
11003013	OPS™ Catheter Pack #13	Sterile	30	1
11003106	OPS™ Catheter Pack #6	Sterile	72	1
11003108	OPS™ Catheter Pack #8	Sterile	72	1

Component details

## Catheter Packs

REORDER	WITH PURIFIL PREFILLED STERILE WATER SYRINGE							
	11003001	11003002	11003009	11003011	11003012	11003013	11003106	11003108
Sterile Field	75cm x 75cm White Plastic	50cm x 50cm Blue Laminated	50cm x 50cm Blue Laminated	70cm x 70cm White Plastic	70cm x 70cm White Plastic	70cm x 70cm White Plastic	60cm x 60cm White Plastic	60cm x 60cm White Plastic
Kidney Dish	2 x 375ml Green					1 x 700ml Clear		
OPS™ Long Tray 24.5cm x 10.5cm x 2.5cm		2 x Teal	2 x Teal	2 x Teal	2 x Teal		1 x Teal	1 x Teal
Gallipots 120ml	2							
Cotton Balls	20	7						5
Gauze Swabs 8ply	2 x 10cm x 10cm		5 x 7.5cm x 7.5cm	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	6 x 10cm x 10cm	6 x 10cm x 10cm
Forceps	4 x Green OPS™	2 x Green OPS™	2 x Green OPS™, 1 x Yellow OPS™	3 x Green	3 x Green			3 x Green
Catheter Securement Device 6.5cm x 10.5cm				1		1		
Lubricating Jelly Sachet 3g				1	1	1	1	1
Sodium Chloride Ampoule				1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml
Prefilled Water Syringe				1 x 10ml	1 x 10ml	1 x 10ml	1 x 10ml	1 x 10ml
Nitrile Powder Free Gloves		2 x Medium					2 x Large	2 x Large
Catheter Insertion Label						14cm x 10cm Yellow		
Paper Towel 40cm x 38cm	2 x 2ply	1 x 2ply	1 x 4ply				1 x 2ply	1 x 2ply
Medical Towel				30cm x 50cm	30cm x 50cm			
Fenestrated Drape 51cm x 67cm (Fen 10cm x 10cm)		1	1					
Fenestrated Drape 60cm x 60cm (Fen 6cm x 10cm)							1 x Blue Cellulose/ Polyester	1 x Blue Cellulose/ Polyester
Fenestrated Drape 80cm x 100cm (Fen 9cm x 18cm)				1	1	1 x SMS		

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

