

Needle Holders - Sterile

Sage Needle Holders are designed to provide surgeons with a secure grip and control, ensuring accurate and efficient suturing. Needle Holders are made of stainless steel, which offers durability and resistance to corrosion. They feature a scissor-like design with a ratcheted locking mechanism that allows for secure needle grasping and holding during suturing.



Finger-ring Handles

Comfortable grip allowing for sustained control



Ratcheting Mechanism

Enables control and secure management of sutures



Blue Handle

Easy to identify as single use



- Suitable for:**
- Hospital and medical
 - Clinics
 - Emergency departments

REF	Size	Feature	Dispenser Box	Sterility	Pcs/Unit	Units/Ctn
11011023	15cm	Mayo Hegar	Yes	Sterile	20	10
11011156	20cm	Mayo	No	Sterile	20	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

