

## Scissors Sharp/Blunt

Forged from stainless steel, Sage Sharp/Blunt Scissors ensure sharpness and resistance to corrosion. These scissors have one sharp tip and one blunt tip for compatibility to different uses.



### Easy Grip

Comfortable shape for a reliable hold even when wearing gloves



### Durable

Stainless steel material provides durability and corrosion resistance



### Ideal for First Aid

Versatility for various medical applications while ensuring patient safety and precision



- Suitable for:**
- First aid
  - Nurse stations
  - Cutting clothing and bandages

REF	Size	Sterility	Pcs/Unit	Units/Ctn
SSS001	125mm	Non-Sterile	50	5

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

PRIMARY BAG



ZIP BAG



BOX

