

Adhesive Aperture Drape

Suresafe Adhesive Aperture Drapes create an absorbent and impervious barrier to fluid strikethrough. They are designed to provide a protective covering to isolate a surgical incision or surgical site. These drapes come in a variety of fenestration sizes and shapes, and with an adhesive strip to secure the drape during the procedure.



Strong and Durable

Made from low linting and abrasion resistant material




Absorbent

Non-woven absorbent outer layer to maintain a dry working area



Impervious

Helps protect against fluid and contaminant strikethrough to aid in infection control



Suitable for:

- Operating theatres
- Emergency departments
- Procedure packs

REF	Size	Hole Size	Hole Shape	Sterility	Pcs/Unit	Units/Ctn
S72261	60cm x 60cm	6cm x 7.5cm	Oval	Sterile	100	1
S72262	75cm x 75cm	6cm x 7.5cm	Oval	Sterile	100	1
S72265	112cm x 150cm	10cm x 15cm	Oval	Sterile	40	1
S72263	90cm x 90cm	10cm x 10cm	Circle	Sterile	100	1
S72281	60cm x 60cm	2.5cm x 6cm	Rectangle	Sterile	200	1
S72282	60cm x 60cm	3.5cm x 10cm	Rectangle	Sterile	160	1
S72283	60cm x 90cm	3.7cm x 11cm	Rectangle	Sterile	160	1
S72284	90cm x 122cm	6cm x 15cm	Rectangle	Sterile	100	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
 1300 995 999 | sales@sentrymedical.com.au
 www.sentrymedical.com.au

