

Fenestrated Drape

Suresafe Fenestrated Drapes create an absorbent and impervious barrier to fluid strikethrough. They are designed to provide a protective covering to isolate a surgical incision or surgical site. These drapes come in a variety of fenestration sizes and shapes.



Strong and Durable

Made from low linting and abrasion resistant material





Absorbent

Non-woven absorbent outer layer to maintain a dry working area



Impervious

Helps protect against fluid and contaminant strikethrough to aid in infection control



Suitable for:

- Operating theatres
- Emergency departments
- Procedure packs

REF	Size	Hole Size	Hole Shape	Sterility	Pcs/Unit	Units/Ctn
S72221	45cm x 60cm	6cm x 7.5cm	Oval	Sterile	120	1
S72222	45cm x 75cm	6cm x 7.5cm	Oval	Sterile	120	1
S72223	60cm x 60cm	6cm x 7.5cm	Oval	Sterile	100	1
S72226	90cm x 90cm	6cm x 7.5cm	Oval	Sterile	100	1
S72230	112cm x 150cm	10cm x 15cm	Oval	Sterile	20	1
S72232	170cm x 270cm	25cm x 25cm	Oval	Sterile	30	1
S72291	60cm x 60cm	2.5cm x 6cm	Rectangle	Sterile	200	1
S72292	60cm x 60cm	3.5cm x 10cm	Rectangle	Sterile	160	1
S72293	60cm x 90cm	3.7cm x 11cm	Rectangle	Sterile	160	1
S72294	90cm x 122 cm	6cm x 15cm	Rectangle	Sterile	100	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

