

Gauze Swabs

Threadx Gauze Swabs are 8-ply cotton swabs ideal for absorbing fluid or being used as either primary or secondary dressings to protect organs from abrasion, drying, or contamination. Threadx Gauze Swabs come French folded to reduce linting and fraying.



Highly Absorbent

Cotton material provides high absorbency




Low Linting

Aids in infection control and minimises the risk of contamination



Flexible

Easy to manipulate and adapt to different shapes and sizes of wounds or surgical sites



GSR2425



GSG2425

Suitable for:

- Absorbing fluids
- Wound padding and protection barrier
- Cleaning and applying topical ointments to wounds

REF	Pcs	Size	Colour	Sterility	Pcs/Unit	Units/Ctn
GSR2223	3pcs	5cm x 5cm	White	Sterile	50	30
GSR2225	5pcs	5cm x 5cm	White	Sterile	50	24
GSR2322	2pcs	7.5cm x 7.5cm	White	Sterile	50	24
GSR2323	3pcs	7.5cm x 7.5cm	White	Sterile	50	24
GSR2325	5pcs	7.5cm x 7.5cm	White	Sterile	50	24
GSG2425	5pcs	10cm x 10cm	Green	Sterile	50	9
GSR2422	2pcs	10cm x 10cm	White	Sterile	50	16
GSR2423	3pcs	10cm x 10cm	White	Sterile	50	8
GSR2425	5pcs	10cm x 10cm	White	Sterile	50	8

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

