

Drainage Bag Accessories Bag Hanger Plastic

Urimaax Drainage Bag Accessories Bag Hanger Plastic is a dual folding plastic hanger ideally used for urinary drainage bags from a catheterised patient.



Durable and Secure

Afixes securely to most bed frame shapes



Rotating Hooks

Round hooks rotate 90 degrees for secure fitting



Easy to Use

Simple and convenient application and removal



- Suitable for:**
- Most bed shapes
 - Urimaax urine drainage bag range

REF	Type	Sterility	Pcs/Unit	Units/Ctn
10013200	Plastic	Non-Sterile	100	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

