

Drainage Bag B/Drain Push-Pull

Urimaax Drainage Bag B/Drain Push-Pull is a large capacity 2000ml drainage bag with a push-pull bottom drain valve and a universal stepped catheter connector. This system has a graduated front with a white frosted backing to allow easy visual readings and hanger holes that fit universal plastic or metal hangers. Ideal for use with Urimaax Drainage Bag Accessories.



Durable and Secure

Kink resistant tubing and strong securement for patient peace of mind



Easy to Read

100ml graduated measurements for easy identification



Fast and Safe

Universal stepped catheter connector and anti-reflux valve



Suitable for:

- Suprapubic or indwelling catheters
- Night catheterisation
- To hold up to 2000ml drainage

| REF | Type | Capacity | Sterility | Pcs/Unit | Units/Ctn |
|----------|--------------|----------|-------------|----------|-----------|
| 10013101 | Push-Pull A1 | 2000ml | Non-Sterile | 200 | 1 |
| 10013102 | Push-Pull A2 | 2000ml | Sterile | 75 | 1 |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

