

Drainage Bag B/Drain T-Bar

Urimaax Drainage Bag B/Drain T-Bar is a large capacity 2000ml drainage bag with a bottom T-bar drain and a universal stepped catheter connector. This system has a graduated front with a white frosted backing to allow easy visual readings and hanger holes that fit universal plastic or metal hangers. Ideal for use with Urimaax Drainage Bag Accessories.



Durable and Secure

Kink resistant tubing and strong securement for patient peace of mind



Easy to Read

100ml graduated measurements for easy identification



Fast and Safe

Universal stepped catheter connector and anti-reflux valve



Suitable for:

- Suprapubic or indwelling catheters
- Night catheterisation
- To hold up to 2000ml drainage

REF	Type	Capacity	Sterility	Pcs/Unit	Units/Ctn
10013103	T-Bar A3	2000ml	Non-Sterile	200	1
10013120	Standard T-Bar A3	2000ml	Non-Sterile	200	1
10013122	T-Bar A2+	2000ml	Sterile	75	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

