

## Drainage Bag Leg

Urimaax Drainage Bag Leg features a fast and safe flip-top bottom drain and a universal stepped catheter connector. The drainage bags have 100ml graduated measurements for easy identification of volume and feature a soft cloth backing for added comfort.



### Durable and Secure

Kink resistant tubing and strong securement for patient peace of mind



### Easy to Read

100ml graduated measurements for easy identification



### Fast and Safe

Flip-top bottom, universal stepped catheter connector, and anti-reflux valve



500ml



750ml

### Suitable for:

- Concealing under clothing
- Enhancing user dignity and comfort
- Use during activities of daily living

REF	Tubing Length	Capacity	Sterility	Pcs/Unit	Units/Ctn
10013302	6cm	500ml	Sterile	80	1
10013303	10cm	500ml	Sterile	80	1
10013304	30cm	500ml	Sterile	80	1
10013305	50cm	500ml	Sterile	80	1
10013402	6cm	750ml	Sterile	80	1
10013403	10cm	750ml	Sterile	80	1
10013404	30cm	750ml	Sterile	80	1
10013405	50cm	750ml	Sterile	80	1
10013503	10m	500ml	Non-Sterile	80	1
10013603	10cm	750ml	Non-Sterile	80	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

