

Drainage Bag Meter 400ml

Urimaax Drainage Bag Meter consists of dual chambers - a 400ml burette and a universal stepped catheter connector. This system has a needle free sampling port and hanger holes that fit universal plastic or metal hangers. Ideal for use with Urimaax Drainage Bag Accessories.



Dual Chambers

For accurate and clearly visible urine measurement



Accessories Included

Dual hangers, hook, bed sheet clamp and rope tie included in pack



Fast and Safe

Universal stepped catheter connector and anti-reflux valve



- Suitable for:**
- Suprapubic or indwelling catheters
 - One-handed application
 - Most bed shapes

REF	Capacity	Sterility	Pcs/Unit	Units/Ctn
10013108	2500ml	Sterile	10	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

