

Drainage Bag Night 2000ml

Urimaax Drainage Bag Night is a large capacity 2000ml drainage bag with an anti-reflux valve and a universal stepped catheter connector. This system has kink resistant tubing with regular bore and hanger holes that fit universal plastic or metal hangers. Ideal for use with Urimaax Drainage Bag Accessories.



Durable and Secure

Kink resistant tubing and strong securement for patient peace of mind



Easy to Read

100ml graduated measurements for easy identification



Fast and Safe

Universal stepped catheter connector and anti-reflux valve



- Suitable for:**
- Night catheterisation
 - Suprapubic or indwelling catheters
 - To hold up to 2000ml drainage

REF	Capacity	Sterility	Pcs/Unit	Units/Ctn
10013110	2000ml	Non-Sterile	250	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

