

Wound Contact Layer

Zorflex® Wound Contact Layer is an activated carbon dressing that provides an alternative to traditional antimicrobial dressings that may contain chemical agents such as silver, iodine or Polyhexamethylene biguanide (PHMB). Naturally antibacterial, antiviral and anti-fungal, Zorflex® is a highly adsorbent, flexible textile which conforms to the body's contours and is easy to cut, fold or layer.



Naturally Antimicrobial

100% pure activated carbon cloth that traps odour and microbes to support wound healing



Conductive Dressing

Restores the body's natural Transepithelial Potential (TEP), aiding healing by encouraging accelerated tissue formation



Protease Modulation

Regulates enzymes in the wound which can promote healing



Moisture Handling

Exudate transfers into the secondary dressing, which can be replaced as necessary



Suitable for:

- Wounds with suspected or proven infection or biofilm presence
- Dry or discharging, partial and full thickness wounds
- Wounds with unviable tissue
- Malodorous wounds

REF	Size	Sterility	Pcs/Unit	Units/Ctn
ZF10	10cm x 10cm	Sterile	10	10
ZF15	15cm x 25cm	Sterile	10	5
ZF20	10cm x 20cm	Sterile	10	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

