


Haemostat Kelly Forceps - Sterile

Forged from stainless steel, Sage Haemostat Kelly Forceps feature horizontal serrations on half of the jaw, a straight end, and a ratchet locking mechanism. These forceps are designed for clamping blood vessels or securely holding heavy tissues.



Ergonomical Handles

Allowing for comfortable and controlled use



Half Serrated Jaws

Providing a strong and reliable grip on the tissue or vessel being clamped



Blue Handle

Easy to identify as single use



- Suitable for:**
- Hospital and medical
 - Clinical
 - Operating suites

REF	Size	Feature	Dispenser Box	Sterility	Pcs/Unit	Units/Ctn
11011028	14cm	Straight	Yes	Sterile	25	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

