


## McIndoe Forceps - Sterile

Forged from stainless steel, Sage McIndoe Forceps feature a slender, elongated design to securely grasp foreign objects during a medical procedure.



### Secure Grip

The textured handles ensure a reliable hold, even when wearing gloves



### Serrated Tips

For safe and secure grasping of tissue or fine objects



### Blue Ends

Serve as a clear indicator that the instrument is a single-use disposable item, ensuring proper identification and adherence to safety protocols



- Suitable for:**
- Hospital and medical
  - Clinics
  - Operating suites

REF	Size	Sterility	Dispenser Box	Pcs/Unit	Units/Ctn
11011021	15cm	Sterile	Yes	30	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

