


## Mosquito Forceps - Sterile

Sage Mosquito Forceps are designed with a slender, elongated shape and fine jaws that taper to a point. The forceps are made of stainless steel, ensuring durability and corrosion resistance.



### Finger-ring Handles

Comfortable grip allowing for sustained control




### Fine Tips

For ease and control when handling small objects or tissues



### Blue Handle

Easy to identify as single use



Straight

Curved



**Suitable for:** - Controlling bleeding in surgical procedures

REF	Size	Feature	Dispenser Box	Sterility	Pcs/Unit	Units/Ctn
11011020	12.5cm	Straight	Yes	Sterile	25	1
11011022	12.5cm	Curved	Yes	Sterile	20	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

