


Semken Tissue Forceps - Sterile

Forged from stainless steel, Sage Semken Tissue Forceps are narrow, toothed forceps designed to handle precise dissection and handling of delicate tissues. They are designed to grasp and hold tissues, sutures, or other objects during surgical operations.



Secure Grip

The textured handles ensure a reliable hold, even when wearing gloves




Extra Fine Serrated Tips

For precise dissection and handling of delicate tissues



Blue Tip

Easy to identify as single use



- Suitable for:**
- Hospital and medical
 - Operating theatres
 - Emergency departments

REF	Size	Sterility	Dispenser Box	Pcs/Unit	Units/Ctn
11011132	12.5cm	Sterile	Yes	50	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

