


Splinter Forceps - Sterile

Forged from stainless steel, Sage Splinter Forceps are designed to securely grasp foreign objects during a medical procedure. These forceps feature fine, pointed and serrated tips to ensure maximum control of the objects.



Secure Grip

The textured handles ensure a reliable hold, even when wearing gloves



Pointed Tips

To handle delicate tissues and foreign bodies



Blue Ends

Serve as a clear indicator that the instrument is a single-use disposable item, ensuring proper identification and adherence to safety protocols



Suitable for:

- Hospital and medical
- Clinics
- Operating suites
- Removal of small foreign objects

REF	Size	Sterility	Dispenser Box	Pcs/Unit	Units/Ctn
11011017	11.5cm	Sterile	Yes	60	10
11011029	12.5 cm	Sterile	Yes	60	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

