

Flaminal® Enzyme Alginogel

Flaminal® is an enzyme alginogel – a wound healing agent that contains a hydrated alginate with an embedded antibacterial enzyme system. It can be used for acute and chronic wound types, or wounds at risk of infection: Flaminal® Hydro for slightly to moderately exuding wounds and Flaminal® Forte for wounds with moderate to high amounts of wound exudate.



Debriding Gel

Dissolves necrotic tissue and fibrinous slough and absorbs them from the wound bed



Absorbent Alginate

Absorbs excess exudate, enables autolytic debridement (cleaning of wound and absorption of bacteria)



Antimicrobial Enzyme System

Offers antimicrobial protection, and reduces bacteria released from a biofilm



Suitable for:

- Superficial and deep partial thickness burns
- Diabetic foot ulcers, leg ulcers or pressure ulcers
- Infected wounds
- Neonates/paediatrics, oncologic, post-surgical or traumatic wounds
- Skin tears

| REF | Description | Sterility | Pcs/Unit | Units/Ctn |
|------|---------------------|-------------|----------|-----------|
| 1015 | Flaminal Forte 50g | Non-Sterile | 1 | 120 |
| 1016 | Flaminal Forte 500g | Non-Sterile | 1 | 18 |
| 1021 | Flaminal Hydro 50g | Non-Sterile | 1 | 120 |
| 1010 | Flaminal Hydro 500g | Non-Sterile | 1 | 18 |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

